



# THE NEW PARO360 SOFTWARE



New, clear paroManager with tablet integration and workflow, new filter options and Multi-database capability.

## 1. MODELLING WITH BOTH FEET

Simple drag and drop of elements from one side to the other side



## 2. CONNECTION TOP AND BOTTOM SIDE PROFILES

Maximum reduction of manual post-processing



## 3. AUTOMATIC GENERATION OF MATCHING 3D BOTTOM PROFILES FROM THE MODELLED SURFACE

All foot shapes can be milled on both sides regardless of shape



## 4. FREE MODELLING OF BOTTOM PROFILES

Individual modelling with all known tools, as well as easy integration of elements, such as reinforcements, cavities etc.



... and many more innovative features.